

Absorbance Reference

Absorbance mode.

OOIBase32 uses an equation to determine the concentration of a species in solution (illustrated below). The software uses this equation to evaluate each pixel on the detector and produce the absorbance spectrum:

$$A_{\lambda} = -\log_{10} \left(\frac{S_{\lambda} - D_{\lambda}}{R_{\lambda} - D_{\lambda}} \right)$$

Where:

S = Sample intensity at wavelength λ

D = Dark intensity at wavelength λ

R = Reference intensity at wavelength λ

The concentration of a species in a solution directly affects the absorbance of the solution. This relationship, known as Beer's Law, is expressed as:

$$A_{\lambda} = \epsilon_{\lambda} c l$$

Where:

A = Absorbance at wavelength λ ,

ϵ_{λ} = Extinction coefficient of the absorbing species at wavelength λ

c = Concentration of the absorbing species and l is the optical path length of the absorption.

Click the Absorbance mode () icon or select **Spectrum | Absorbance Mode** to enter Absorbance mode.

As a reference, here some calculated numbers at 266nm

Protein	Space group	conc(M)	eps (M⁻¹cm⁻¹)	Absorbance 100µm crystal
elastase	P212121	0.0309	45410	14
insulin	I213	0.0843	4820	4
hewl	P43212	0.0610	32717	20
pyp	P65	0.0592	10518	6
rnaseA	P3221	0.0433	6998	3
thaumatin	P41212	0.0592	24823	15
trypsin	P212121	0.0315	30478	10

Abs=conc*eps*l , l=path length crystal