## Coherent Diffraction Imaging at Third and Fourth Generation X-Ray Sources

<u>Chapman H. N.</u><sup>1</sup>, Barty A.<sup>1</sup>, Beetz T.<sup>5</sup>, Cui C.<sup>2</sup>, Hajdu J.<sup>6</sup>, Hau-Riege S.P.<sup>1</sup>, He H.<sup>3</sup>, Howells M.<sup>2</sup>, London R.<sup>1</sup>, Marchesini S.<sup>1</sup>, Shapiro D.<sup>4</sup>, Spence J.C.H.<sup>3</sup>, Weierstall U.<sup>3</sup>

<sup>1</sup>University of California, Lawrence Livermore National Laboratory, Livermore, CA 94550, USA, Email: henry.chapman llnl.gov, <sup>2</sup>Advanced Light Source, Lawrence Berkeley National Laboratory, Berkeley, CA 94720, USA, <sup>3</sup>Department of Physics and Astronomy, Arizona State University, Tempe, AZ 85287, USA, <sup>4</sup>Center for Biophotonics Science and Technology, UC Davis, Davis CA 95817, USA, <sup>5</sup>Department of Physics and Astronomy, Stony Brook University, Stony Brook, NY 11794, USA, <sup>6</sup>Biomedical Center, Uppsala University, Box 596, SE 751 24, Uppsala, Sweden

We are assessing X-ray diffraction microscopy as a means to perform high-resolution three-dimensional characterisation of non-periodic isolated objects (particles). Several recent experimental and computational developments have enabled us to perform full 3D X-ray diffraction imaging, with high resolution in all three dimensions. These 3D reconstructions were performed from the diffraction data alone. These developments include the Stony Brook diffraction apparatus [1], which allows 3D diffraction datasets to be quickly acquired; the Shrinkwrap phase-retrieval algorithm [2], which allows images to be reconstructed *ab initio* from incomplete diffraction datasets; and a fast distributed FFT [3] and reconstruction software implemented on a computer cluster, which allows  $1024^3$ diffraction datasets to be phased in several hours. We have achieved high-resolution 3D reconstructions of both well-characterized test objects and of mesoporous foam particles that cannot be otherwise characterized. We find that high resolution imaging of thick objects can only be attained in the context of 3D measurement and reconstruction. Reconstruction from diffraction data acquired over many sample orientations allows one to avoid defocus (depth of field) artifacts as well perform a quantitative measurement of refractive index that is not possible from single-view diffraction data.

Resolution of X-ray diffraction imaging will ultimately be limited by radiation damage [4]. One eventual goal is to surpass damage resolution limits of individual particles by using streams of identical particles, such as protein macromolecules, using flash-imaging by X-ray free-electron lasers (XFELs) [5]. Models show that these methods should allow close to atomic resolution imaging. Upcoming experiments at the DESY VUV-FEL will be described where we will attempt to perform diffraction imaging of the damage of particles as they are exposed to short, intense X-ray pulses, to verify models of short-pulse damage.

This work was performed under the auspices of the U. S. DOE by LLNL under Contract No. W-7405-ENG-48 and by LBNL under Contract No. DE-AC03-76SF00098. The Center for Biophotonics, an NSF Science and Technology Center, is managed by UC Davis, under Contract No. PHY0120999. Additional support was provided by the Swedist Research Council.

## References

- T. Beetz, M. R. Howells, C. Jacobsen, C.-C. Kao, J. Kirz, E. Lima, T. O. Mentes, H. Miao, C. Sanchez-Hanke, D. Sayre, D. Shapiro, Nucl. Instrum. Meth Phys. Res. A (2005, in press).
- [2] S. Marchesini, H. He, H. N. Chapman, S. P. Hau-Riege, A. Noy, M. R. Howells, U. Weierstall, and J. C. H. Spence, Phys. Rev. B 68, 140101(R) (2003).
- [3] R. Crandall et al. http://images.apple.com/acg/pdf/20040827\_GigaFFT.pdf (2004).
- [4] M.R. Howells et al, J. Electron. Spectrosc. Rel. Phenom (2005, submitted).
- [5] R. Neutze, R. Wouts, D. van der Spoel, E. Weckert, and J. Hajdu, Nature 406 753-757 (2000).