

Toward *In Situ* Mapping of Crystal Selection during OFZ Crystal Growth

Peter Khalifah^{1,2}, Jonathan Denney¹, Yusu Wang¹, Adam Corrao¹

Simon Billinge^{2,3}, Christopher J. Wright³, Songsheng Tao³, Chia-Hao Liu³

Eric Dooryhee², David Sprouster^{1,2}

W. Adam Phelan⁴, Lucas A. Pressley⁴

Katsuyo Thornton⁵, Guanglong Huang⁵, Mojue Zhang⁵, Praveen Soundararajan⁵, David Taboada⁵

¹Stony Brook University, ²Brookhaven National Laboratory, ³Columbia University,

⁴Johns Hopkins University (PARADIM), ⁵University of Michigan

This work was supported as part of GENESIS: A Next Generation Synthesis Center, an Energy Frontier Research Center funded by the U.S. Department of Energy, Office of Science, Basic Energy Sciences under Award Number DE-SC0019212



U.S. DEPARTMENT OF
ENERGY

Office of
Science



BROOKHAVEN
NATIONAL LABORATORY

OAK RIDGE
National Laboratory

Colorado
State
University



UC San Diego

COLUMBIA UNIVERSITY
IN THE CITY OF NEW YORK


Farmingdale
State College
State University of New York

M UNIVERSITY OF
MICHIGAN



A NEXT **GENERATION**
SYNTHESIS CENTER

Large single crystals have been transformative



- Source of many key insights into high- T_c cuprate superconductors (inelastic neutron, ARPES, etc.)
- More recently being applied to other fields (solar fuels, catalysis, etc.)
- Key technology for producing crystals is optical floating zone furnace
- **OFZ growth is currently an art and not a science** – a quantitative understanding is lacking!

Some specific GENESIS goals

- Apply modern synchrotron methods to probe synthesis *in situ* in new ways and with unprecedented sensitivity
- Simultaneously probe polycrystalline, liquid, and crystal regions
- Obtain new insights into and greater control over crystal growth processes
- Gain access to important hard-to-grow phases



U.S. DEPARTMENT OF
ENERGY

Office of
Science



BROOKHAVEN
NATIONAL LABORATORY

OAK RIDGE
National Laboratory

Colorado
State
University



UC San Diego

COLUMBIA UNIVERSITY
IN THE CITY OF NEW YORK

Farmingdale
State College
State University of New York

M UNIVERSITY OF
MICHIGAN



A NEXT **GENERATION**
SYNTHESIS CENTER

Large single crystals have been transformative



- Source of many key insights into high- T_c cuprate superconductors (inelastic neutron, ARPES, etc.)
- More recently being applied to other fields (solar fuels, catalysis, etc.)
- Key technology for producing crystals is optical floating zone furnace
- **OFZ growth is currently an art and not a science** – a quantitative understanding is lacking!

Some specific GENESIS goals

- Apply modern synchrotron methods to probe synthesis *in situ* in new ways and with unprecedented sensitivity
- Simultaneously probe polycrystalline, liquid, and crystal regions
- Obtain new insights into and greater control over crystal growth processes
- Gain access to important hard-to-grow phases



U.S. DEPARTMENT OF
ENERGY

Office of
Science



BROOKHAVEN
NATIONAL LABORATORY

 **COLUMBIA UNIVERSITY**
IN THE CITY OF NEW YORK

 **OAK RIDGE**
National Laboratory

 **Farmingdale State College**
State University of New York


 **Colorado State University**
 **UNIVERSITY OF MICHIGAN**



UC San Diego
 **BERKELEY LAB**

A NEXT **GENERATION**
SYNTHESIS CENTER

Opportunities to probe OFZ crystal growth



Sample size:

- ~5 mm diameter
- ~5 cm length

Crystal growth in OFZ furnace typically monitored using only video

- Can roughly distinguish between solid and liquid regions
- Cannot monitor temperature / temperature profiles
- Cannot probe the interior of the crystal growth rods

New opportunities using a synchrotron, but **boules are not standard samples!**



U.S. DEPARTMENT OF
ENERGY

Office of
Science



BROOKHAVEN
NATIONAL LABORATORY

OAK RIDGE
National Laboratory

Colorado
State
University



UC San Diego

COLUMBIA UNIVERSITY
IN THE CITY OF NEW YORK

Farmingdale
State College
State University of New York


M UNIVERSITY OF
MICHIGAN



A NEXT **GENERATION**
SYNTHESIS CENTER

(Mis)utilization of single crystal diffraction

~0.5 mm crystal



~10 mm crystal



Routine to characterize small xtals (< 0.5 mm) with SXRD to get 3 x 3 UB matrix

Very hard to characterize large crystals (5 – 10 mm) by Laue or SXRD

- Laue – works poorly when symmetry is low
- SXRD – transmission of beam is essentially zero; crystal larger than beam

(Mis)utilized high-energy synchrotron X-rays



Experiments done at 28-ID powder diffraction beamlines of NSLS-II

- Suitable for mapping studies - detector / motors limit speed (≤ 1 sec)
- High energy X-rays can penetrate through (some) 5 mm thick samples
- Ongoing construction of HEX beamline (photons to ~ 200 keV)



U.S. DEPARTMENT OF
ENERGY

Office of
Science



BROOKHAVEN
NATIONAL LABORATORY

OAK RIDGE
National Laboratory

Colorado State
University



UC San Diego

COLUMBIA UNIVERSITY
IN THE CITY OF NEW YORK

Farmingdale State College
State University of New York

M UNIVERSITY OF MICHIGAN



A NEXT **GENERATION**
SYNTHESIS CENTER

1. Investigating the crystal selection process

Ideal growth product is a single crystal


- What is the pathway from powder to single crystal?
- What growth variables help/hinder selection process?
- Knowledge should improve crystal size and quality

It is possible but destructive / laborious to study *ex situ*

- What are opportunities to study this process *in situ*?
- Can optimize growths more rapidly when using *in situ* data

Initial *ex situ* synchrotron studies have been carried out

- Will be working to extend methods to *in situ* experiments
- May be applicable to processes other than just OFZ crystal growth



U.S. DEPARTMENT OF
ENERGY

Office of
Science



BROOKHAVEN
NATIONAL LABORATORY

OAK RIDGE
National Laboratory

Colorado State
University



UC San Diego

COLUMBIA UNIVERSITY
IN THE CITY OF NEW YORK

Farmingdale State College
State University of New York

UNIVERSITY OF MICHIGAN



A NEXT **GENERATION**
SYNTHESIS CENTER


Ex situ data collection on powder beamline



Transmission measurements on TiO_2 boule grown in OFZ furnace

- Utilized high-energy X-rays (67 keV) and a 2D detector
- Sample scanned using a medium grid (0.2 x 0.2 mm spot size & step size)
 - Probe region about 5 mm x 40 mm, so about $25 \times 400 = 10,000$ grid points
- No rotation of sample was done in these proof-of-principle experiments
 - OFZ furnaces are designed to rotate samples during growth


Variation in diffraction patterns from crystal boule



Diffraction patterns evolve from those of single crystal to multi-crystal to powder

- Existing software designed to deal with crystal spots or powder rings but not both

Analysis of crystal growth product - grain tracking



Single crystal diffraction spots used to extract info about grains in boule

- When beam scanned within a single crystal domain, same diffraction spots are seen
- For any one hkl reflection at a given azimuthal angle, can make a map of the grain
- Utilized existing peak finding functionality of trackpy package

Simple but effective approach

Analysis of crystal growth product - grain tracking





Image shows crystal domains which produce the 10 most intense diffraction spots

- Can see single crystal domains both expanding and terminating
- All major grains appear of have nucleated near the center of the boule
- Unexpected intensity bands commonly seen – ascribed to mosaicity effects

Next step: Calculate orientation-dependence of TiO_2 growth rate

Rocking curve measured at single spot




Rotated crystal by 1° in 0.01° steps

- Sharp onset and loss of intensity even with 0.01° rotation
- Presently testing precession methods to avoid intensity bands
- Suggests utility for investigating influence of growth parameters on quality

Mosaic spread variation across crystal



Boule 1



Can also probe mosaic spread in context of maps

- Can extract position-dependent variation in azimuthal angle (right)
- More data needed to extract full mosaicity matrix (more spots, +rotation)
- Working to develop infrastructure needed to study crystal growth *in situ*



U.S. DEPARTMENT OF
ENERGY

Office of
Science



BROOKHAVEN
NATIONAL LABORATORY

OAK RIDGE
National Laboratory

Colorado
State
University

UC San Diego

A NEXT **GENERATION**
SYNTHESIS CENTER


COLUMBIA UNIVERSITY
IN THE CITY OF NEW YORK

Farmingdale
State College
State University of New York


M UNIVERSITY OF
MICHIGAN

BERKELEY LAB


Development of infrastructure for *in situ* studies



Mark I – air only



Mark II – vacuum / gas





Mark III – dual-rod

Carrying out design and fabrication of “mini-FZ” furnaces for *in situ* studies

- Avoids limitations of commercial FZ furnaces (high cost, large mass)

Development of infrastructure for *in situ* studies




Have used powder diffraction data to extract temperature profiles of heated rods

- From thermal expansion behavior (strain), can get T with precision < 0.25 °C
- Have developed physical models (dashed lines) that reproduce results
- First model for FZ environment parameterized and validated from expt. data

Next steps: (1) Inline analysis, (2) *in situ* crystal growth, (3) automation of growth

Automated data acquisition at NSLS-II, BNL



U.S. DEPARTMENT OF
ENERGY

Office of
Science



BROOKHAVEN
NATIONAL LABORATORY

OAK RIDGE
National Laboratory

Colorado
State
University

UC San Diego

COLUMBIA UNIVERSITY
IN THE CITY OF NEW YORK





Farmingdale
State College
State University of New York

M UNIVERSITY OF
MICHIGAN



A NEXT GENERATION
SYNTHESIS CENTER

Analysis server for automated grain tracking



Conclusions / Acknowledgements

- Actively developing synchrotron methods to probe crystal growth
 - Temperature distributions <https://doi.org/10.1107/S1600576720007062>
 - Crystal grain distributions <https://pubs.acs.org/doi/pdf/10.1021/acs.chemmater.1c00602>
- Methods are or will be compatible with *in situ* crystal growth

DOE support: This work was supported as part of GENESIS: A Next Generation Synthesis Center, an Energy Frontier Research Center funded by the U.S. Department of Energy, Office of Science, Basic Energy Sciences under Award Number DE-SC0019212

NSF support: This material is based upon work supported by the National Science Foundation (Platform for the Accelerated Realization, Analysis, and Discovery of Interface Materials (PARADIM)) under Cooperative Agreement No. DMR-1539918.

User facility support: This research used the National Synchrotron Light Source-II, a U.S. Department of Energy (DOE) Office of Science User Facility operated for the DOE Office of Science by Brookhaven National Laboratory under Contract No. DE-SC0012704.